



Optistat™ Handheld CT Contrast Delivery System

Covidien continues to reinforce its position as a total solutions provider for contrast delivery with the Optistat™ CT handheld contrast delivery system. The first ultra compact system of its kind, it allows single-handed power injection with excellent economy, space efficiency, portability and control. It is ideal for CT, x-ray and mobile coaches. It can also be used as a backup injector where needed.

Order Information

| Description | Qty | Unit | Size | Order # | 11 Digit NDC |
|---|-----|------|-----------------------|---------|--------------|
| Console for handheld CT contrast delivery system | 1 | each | | 810025 | |
| CT contrast delivery system, handheld | 1 | each | | 810001 | |
| CT contrast delivery system, handheld, with console option | 1 | each | | 810000 | |
| Interface kit, Optistat™ handheld contrast delivery system to CT 9000™ ADV system | 1 | each | | 810620 | |
| Interface kit, Optistat™ handheld contrast delivery system to Optistat™ system | 1 | each | | 810615 | |
| Mount, IV pole for CT handheld powerhead | 1 | each | | 810190 | |
| Mount, rolling stand for CT handheld powerhead | 1 | each | | 810175 | |
| Mount, wall for CT handheld powerhead | 1 | each | | 810185 | |
| Y-tubing adaptor for connecting Optistat™ contrast delivery system to CT 9000™ ADV contrast delivery system | 50 | each | 6" (each adaptor leg) | 810555 | |

Features and Benefits

| Feature | Benefit |
|---|--|
| Accommodates Ultraject™ 125 mL prefilled syringes | Use of Ultraject™ prefilled syringes help reduce the potential for medication errors ¹ , help meet stringent regulatory requirements ¹ , may reduce the risk of infection ¹ , and improve efficiency ² . |
| Optional interface with another Optistat™ handheld contrast delivery system or a CT 9000™ ADV contrast delivery system | Excellent economical option for dual head protocols. |
| Small and lightweight | Easily portable for multiple uses within the radiology department, eliminating hand injections for IVPs, pediatrics, etc. |
| Low cost alternative | Ideal for emergency backup situations |

Easy to install

Installs in most applications in less than an hour with minimal space requirements

Frequently Asked Questions

Question

Answer

What syringe sizes can be used with the Optistat™ handheld contrast delivery system?

125 mL Ultraject™ prefilled syringes and 130 mL empty CT syringes when manually filled

What mounting options are available for the Optistat™ handheld contrast delivery system?

IV pole and wall hanger mounts provide flexible options

Instructions for use may be found in the Optistat™ Operator's Manual.
